## PRÁCTICA I.2

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## EXERCISE I.2

### DETERMINATION OF LIPOPHILICITY SUBSTITUENT CONSTANTS OF SULFONAMIDES BY MEANS OF REVERSED- PHASE THIN-LAYER CHROMATOGRAPHY

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*E-mail: doraprad@dedalus.lcc.ufmg.br*

#### INTRODUCTION

The partition coefficient ( $P$ ) is one of the most important factors in controlling drug action in biological systems. Hansch and coworkers [1] derived substituent constants for the contribution of individual atoms or groups to the partition coefficient, the lipophilicity (or hydrophobicity, or Hansch) constant ( $\pi$ ). This is defined as the  $\log P_{RX} - \log P_{RH}$  (eq. 1), where  $P_{RX}$  and  $P_{RH}$  are the partition coefficients, determined in system 1-octanol-water, of substituted and unsubstituted compounds, respectively. Thus, the constant  $\pi$  indicates the change in the logarithm of the partition coefficient resulting from the introduction of a substituent  $X$ .

$$\pi_x = \log \frac{P_{RX}}{P_{RH}} = \log P_{RX} - \log P_{RH}$$

where  $P_{RX}$  = partition coefficient of  $RX$ ,  $P_{RH}$  = partition coefficient of  $RH$ .

The  $\pi$  values for many substituent groups were directly measured for a variety of drugs and used in order to calculate the partition coefficient of other compounds. However, the direct determination of partition coefficient is tedious and often presents practical difficulties, particularly when the compound is highly insoluble in either of the solvent phases. In order to avoid these practical difficulties, Boyce and Milborrow [2] had proposed the chromatographic  $Rm$  value as an expression of lipophilic character of molecules. The  $Rm$  is calculated using the  $Rf$  value determined by means of reversed-phase thin-layer chromatography. Equation 2 is used to calculate the  $Rm$  value for neutral substances and eq. 3 for acids and bases [3]. The change in the value of  $Rm$  ( $Rm_{RX} - Rm_{RH}$ ) for a substituent is a free-energy-based constant identical with  $\pi$  used by

Hansch (eq. 4).

$$Rm = \log\left(\frac{1}{Rf} - 1\right) \quad (2)$$


$$Rm = \log\left(\frac{1}{Rf} - 1\right) + \log\left(\frac{Ka + [H^+]}{[H^+]}\right) \quad (3)$$

where  $Ka$  = dissociation constant,  $[H^+]$  = hydronium ion concentration of the mobile phase

$$\pi_x = Rm_{RX} - Rm_{RH} \quad (4)$$


The simple laboratory activity described here (reverse-phase thin-layer chromatography of antibacterial drugs, the sulfonamides) offers an opportunity to show the differences between the normal and reverse-phase thin-layer chromatography and introduces a new chromatographic parameter, the  $Rm$ . By means of this experiment, it is possible to determine  $\pi$  values of three substituent groups present in the structures of sulfonamides and calculate the partition coefficient of the compounds.

**R= H sulfanilamide pKa= 10,45**

**R=  sulfamethoxypyridazine pKa= 7,05**

**R=  sulfathiazole pKa= 7,10**

**R=  sulfamethazine pKa= 7,70**

## EXPERIMENTAL

Basic references: G. L. Biagi, A. M. Barbaro, M. C. Guerra, G. C. Forti, M. E. Fraccasso. *J. Med. Chem.* **17**, 28 (1974); M. A. F. Prado. *J. Chem. Educ.* **78**, 533 (2001).

- Stationary phase: Mix 6 g of silica gel with 15 mL of water, shake, coat the plates (12 glass plates measuring 9 × 4.5 cm; 0.25 mm), dry at room temperature, heat in oven at 120 °C for 4 h, and cool at room temperature. A stationary nonpolar phase is obtained by impregnation of the silica gel layer with a 5 % 1-octanol in ethyl ether. The impregnation is carried out by immersion of the plates in the solution. The silica

- gel layer is removed at 3 mm from each lateral edge and 1 cm from the high edge.
- Mobile phase: An aqueous buffer at pH 7.4 (mixture of 50 mL of 0.1 mol/L potassium dihydrogen phosphate solution and 19.5 mL of 0.1 mol/L sodium hydroxide solution).
  - Samples: The sulfonamides (sulfanilamide, sulfathiazole, sulfamethazine, and sulfamethoxypyridazine) are dissolved in acetone (10 mg in 25 mL).
  - Chromatography: The solutions of sulfonamides are spotted on a line 1 cm from the lower edge of the plate using a glass capillary. The plate is left in the chromatographic chamber until the mobile phase reaches the end of the stationary layer.
  - Revelation: The developed plate is dried with hot air, and the spots appear by spraying the plate with a solution prepared by adding 1 mL of concentrated hydrochloridric acid to 100 mL of 0.1 % *p*-dimethylaminebenzaldehyde solution in ethanol.

### CALCULATIONS

- Calculate the *R<sub>f</sub>* value of each sulfonamide.
- Calculate the *R<sub>m</sub>* value of each sulfonamide.
- Calculate the lipophilicity substituent constant ( $\pi$ ) of substituent groups of sulfathiazole, sulfamethazine, and sulfamethoxypyridazine.
- The partition coefficient of sulfanilamide is 0.15. Calculate the partition coefficient of sulfathiazole, sulfamethazine, and sulfamethoxypyridazine.

### REFERENCES

1. P. Hansch, P. P. Maloney, T. Fujita. *Nature* **194**, 178 (1962).
2. C. B. C. Boyce, B. V. Milborrow. *Nature* **208**, 537 (1965).
3. G. L. Biagi, A. M. Barbaro, M. C. Guerra, G. C. Forti, M. E. Fraccasso. *J. Med. Chem.* **17**, 28 (1974).

**Maria Auxiliadôra Fontes Prado**

doraprad@dedalus.lcc.ufmg.br

High standards in safety measures should be maintained in all work carried out in Medicinal Chemistry Laboratories. The handling of electrical instruments, heating elements, glass materials, dissolvents and other inflammable materials does not present a problem if the supervisor's instructions are carefully followed.

This document has been supervised by Prof. Maria Auxiliadôra Fontes Prado) who has informed that no special risk (regarding toxicity, inflammability, explosions),

outside of the standard risks pertaining to a Medicinal Chemistry laboratory exist when performing this exercise.

If your exercise involves any "special" risks, please inform the editor.