PRÁCTICA I.6

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
E-mail: rcattana@exa.unrc.edu.ar

Absorption, distribution, biotransformation, and excretion of a medicinal drug imply its penetration of a cellular membrane. Thus, it is essential to know the mechanism by which drugs pass through the membranes as well as the physicochemical properties of these molecules and the membranes that affect such transference.


The size and molecular shape, the solubility at the absorption site, the degree of ionization, and the relative solubility in lipids of the ionized and nonionized forms are important properties to determine in these drugs.

Several physicochemical properties of a series of compound derivatives of 6,7-diaryl pteridines (Scheme 1), which have been shown to have in vitro antihelmintic properties [1], have been studied.


- 1: R₁ = methyl, R₂ = *p*-Br-phenyl X = Y = O
- 2: R₁ = methyl, R₂ = *p*-F-phenyl, X = Y = O
- 3: R₁ = methyl, R₂ = *p*-Cl-phenyl, X = Y = O
- 4: R₁ = H, R₂ = *p*-methoxyphenyl, X = O, Y = S
- 5: R₁ = H, R₂ = *p*-F-phenyl, X = O, Y = S
- 6: R₁ = H, R₂ = *p*-Cl-phenyl, X = O, Y = S


7: R₁ = H, R₂ = *p*-Br-phenyl, X = Y = S

8: R = 2-thienyl, X = S, Y = H

9: R = *p*-F-phenyl

10: R = *p*-Cl-phenyl

11: R = *p*-Br-phenyl

12: R = 2-thienyl, X = Y = NH₂


13: R = 2-thienyl, X = NH₂, Y = H

14: R = phenyl, X = Y = NH₂

15: R = phenyl, X = NH₂, Y = H

Scheme 1

Some of the studied compounds can show functional tautomerism in which the NH₂ group is involved, such as is shown in Scheme 2.

Scheme 2

On the other hand, it is also possible that these compounds in solution form aggregates through hydrogen bonding.

In order to relate the therapeutic activity of these compounds with the experimentally measured physicochemical properties, it is necessary to know which is the predominant tautomeric form and which are free or aggregated in solution.

Consequently, studies were performed to determine the effects of temperature and concentration in the UV/vis spectra [2] in the solvents and concentrations similar to the ones used for the determination of its bioactivity.

Problem 1

Describe what you expect to observe in the UV/vis spectra by changing temperature and/or concentration in the case of (a) aggregates formation and (b) different forms in equilibrium due to intramolecular formation of tautomers.

Many studies have shown that the biological activity of the substances is predominantly a function of their hydrophobic nature, which is the tendency to repel water and be miscible in nonpolar solvents (lipophilicity).

The hydrophobic/lipophilic nature of a compound is a complex physicochemical property that can be determined by dissolving it in a two-phase mixture of water and an immiscible liquid and then measuring its equilibrium concentration in each phase. Thus, the partition coefficient of the compound between water and the solvent is calculated

The most useful parameter as a measure of hydrophobicity is the partition coefficient between water and octanol ($P_{o/w}$) or to use it in free energy lineal relationships as its logarithm ($\log P_{o/w}$) [3].

Alternatives, particularly chromatographic methods, have been evaluated and used successfully to assess the lipophilicity of organic compounds. Reverse phase liquid chromatography (RPLC) is particularly useful for estimating the physicochemical properties of organic compounds and is known as “structure–retention quantitative relationships”, which are also related to “structure–activity quantitative relationships”. The lipophilicity parameter is obtained experimentally as the logarithm of the capacity factor ($\log k'$)[†] which is proportional to the logarithm of partition coefficient. The use of RPLC to determine lipophilicity descriptor has several advantages over the traditional methods of obtaining partition coefficients because it is faster, no purification of the sample is needed, and it can be performed automatically. One of the most used approximations to determine capacity factors as lipophilicity descriptors is the $\log k'_w$ (logarithm of the capacity factor in pure water). The advantage of using this descriptor is the fact that it is independent of the organic solvent used in its determination, and reflects the partition between the polar and nonpolar media in a way similar to a partition coefficient but is independent of the solute structure.

The $\log k'_w$ can be estimated by the intercept of the plot of eq. 1 [4]

$$\text{Log } k' = S \% + \log k'_w \quad (1)$$

where S is the Snyder force and % is the percentage of the mobile phase.

Equation 2 can also be calculated

$$\text{Log } k' = a + b \cdot E_T \quad (30) \quad (2)$$

[†]Capacity factor = $k' = \frac{\text{Time of peak maximum} - \text{Time of baseline}}{\text{Time of baseline}}$, measured in minutes of the solvent.

where k' is the capacity factor of a given composition in the mobile phase and $E_T(30)$ is the Reichardt's solvent polarity parameter [5] while a and b are the regression coefficients. The value of $\log k'_w$ in pure water can be calculated as the value of $E_T(30)$ corresponding to pure water. It has been demonstrated [2,4] that the values of $\log k'_w$ obtained by means of eq. 2 are more exact than the ones obtained by plotting $\log k'$ vs. % of the organic solvent (eq. 1) and extrapolating to zero concentration of this solvent.

Problem 1

Table 1 shows the values of $\log k'$ measured in MAI columns (artificially immobilized membranes) in mixtures of water-acetonitrile at physiological pH as mobile phase, % and $E_T(30)$ for different solvent compositions of compounds 4, 9, and 15.

Using any statistics software calculates $\log k'_w$ using the described methods (equations 1 and 2). Compare the validity of the values obtained by both methods.

Problem 3

Table 2 shows the values of $\log k'_w$ obtained for all the studied compounds Error! Bookmark not defined. (except those that should be calculated in problem 2) together with the values of $-\log IC_{50}$ Error! Bookmark not defined.

Relate the values of $(-\log IC_{50})$ with those of $(\log k'_w)$. What conclusions can be drawn regarding the structure-activity relationship?

Table 1 Values of % of acetonitrile (ACN), E_T (30), and log k'.

% ACN	Log k' (4)	Log k' (9)	Log k' (15)	E _T (30)
0				63.4
20	1.19	1.57	1.30	59.60
30	0.32	0.81	0.39	58.36
35	-0.017	0.48	0.057	57.74
40	-0.37	0.15	-0.35	57.29
45	-0.68	-0.04	-0.71	56.87

Table 2 Values of antihelminthic activity (-log IC₅₀, ref. 1)^a and logarithm of the capacity factor in pure water in columns MAI (log k'_w) for the compounds in Scheme 1.

Compound	-Log IC ₅₀	Log k' _w
1	2.52	4.84
2	2.10	3.63
3	3.00	5.04
4	1.19	
5	0.72	3.68
6	2.12	4.72
7	3.00	5.68
8	1.08	3.60
9	1.52	
10	2.46	5.26
11	2.71	5.25
12	1.11	3.42
13	0.60	3.62
14	1.22	3.85
15	0.97	

^aIC₅₀ is the inhibitory concentration required to obtain up 50 % in the reduction of population growth of nematode.

REFERENCES

1. C. Ochoa, J. Rodríguez, M. L. Lopez García, A. R. Martínez, M. M. Martínez. Arzneim. *Forsch. Drug. Res.* **46**, 643 (1996).
2. M. Reta, L. Giacomelli, M. Santo, R. Cattana, J. J. Silber, C. Ochoa, M. Rodríguez, A. Chana. *Biomed. Chromatogr.* **17**, 365 (2003).
3. R. Kalizan, A. Nasal, A. Bucinski. *Chem. Rev.* 91 (1993).
4. J. J. Michels, J. G. Dorsey. *J. Chromatogr.* **499**, 435 (1990)
5. C. Reichardt. *Solvents and Solvent Effects in Organic Chemistry*, 2nd ed., VCH, New York (1990).

Rosa Cattana

rcattana@exa.unrc.edu.ar

High standards in safety measures should be maintained in all work carried out in Medicinal Chemistry Laboratories. The handling of electrical instruments, heating elements, glass materials, solvents and other inflammable materials does not present a problem if the supervisor's instructions are carefully followed.

This document has been supervised by Prof. Dr. Rosa Cattana (rcattana@exa.unrc.edu.ar) who has informed that no special risk (regarding toxicity, inflammability, explosions), outside of the standard risks pertaining to a Medicinal Chemistry laboratory exist when performing this exercise.

If your exercise involves any "special" risks, please inform the editor.