PRÁCTICA I.8

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E-mail: madalena@ff.up.pt

INTRODUCTION

Molecular recognition plays an important role in several vital processes, such as protein synthesis, enzymatic catalysis, DNA replication, viral replication, and drug action. On the other hand, artificial molecular recognition can result from a mimetization of the necessary conditions for molecular recognition by nonbiological compounds with structure and functions similar to the biological counterparts [1]. The formation of a complex of the “small molecule-macromolecule”, involving several interactions (hydrophobic, van der Waals, dipole-dipole, hydrogen bonding, charge transfer, etc.), is necessary for the occurrence of a biological function. The formation of such a complex also occurs in artificial recognition, and if there is an entrance of the small molecule into the cavity of the nonbiological macromolecule the new chemical entity is called “inclusion complex”. Examples of artificial molecular recognition include the formation of inclusion complexes between drugs and cyclodextrins, naturally occurring macrocyclic compounds composed of α -D-(+)-glucopyranose units, obtained by enzymatic degradation of starch. There are three particularly important cyclodextrins: α -cyclodextrin, β -cyclodextrin, and γ -cyclodextrin, containing 6, 7, and 8 α -D-(+)-glucopyranose units, respectively. The shape of these compounds is very peculiar, frequently called a “basket” or a “doughnut” with one rim larger than the other. One rim of the cavity possesses secondary hydroxyl groups, whilst the other rim is

functionalized with primary hydroxyl groups. The interior of the cavity is nonpolar, and the external surface and the rims are polar. The interaction of the small molecules with the cavity of cyclodextrins is normally hydrophobic. However, the substituents of the small molecules can also form a hydrogen bonding with the rims of cyclodextrin molecules [2].

Inclusion complexes can be made in an aqueous solution or in a solid phase, and can be analyzed by spectrophotometric methods (IR, UV), nuclear magnetic resonance spectroscopy (^1H and ^{13}C), X-ray diffractometry, and circular dichroism [3–8].

AIM

Preparation of the inclusion complexes of diazepam using α -, β -, and γ -cyclodextrins by solid phase.

MATERIAL AND METHODS

Materials

Diazepam was obtained from Sigma Chemical Co., γ -cyclodextrin was purchased from Fluka, α - and β -cyclodextrins were obtained from Tokyo TCI Kasei.

IR spectra were recorded in ATI Mattson Genesis series FTIR, in 1 % potassium bromide disks (Uvasol, Merck).

Methods [3]

Weigh 5.7 mg of diazepam, 29.2 mg of α -cyclodextrin, 34.0 mg of β -cyclodextrin, and 38.9 mg of γ -cyclodextrin.

Physical mixtures: mix appropriate amount of each cyclodextrin with diazepam. Record the IR spectra.

Inclusion complexes: Grind the remaining amount of physical mixture in quartz mortar during 10 min and then record the IR spectra of the complexes.

RESULTS AND DISCUSSION

Inclusion complexes of diazepam are formed by the solid-phase method. IR spectrophotometry is used to determine which cyclodextrin inclusion complex is formed. Figure 1 shows the IR spectra of the physical mixtures and the inclusion products of diazepam with α -, β -, and γ -cyclodextrins in the carbonyl stretching region.

In the case of γ -cyclodextrin inclusion complex, the 1685 cm^{-1} band was shifted to 1665 cm^{-1} , suggesting the formation of the intermolecular hydrogen bonding between diazepam and γ -cyclodextrin. The IR spectra showed that the inclusion complexes with the diazepam and γ -cyclodextrins were not formed.


Fig. 1 IR spectra of diazepam- α -cyclodextrin (**A**), diazepam- β -cyclodextrin (**B**), and diazepam- γ -cyclodextrin (**C**) physical mixtures (—) and inclusion complexes (---)

REFERENCES

1. P. D. Beer, P. A. Gale, D. K. Smith. *Supramolecular Chemistry*, pp. 57 and 92, Oxford University Press, New York (1999).
2. <<http://www.betadexcyclodextrin.com/html/properties.html> (07/01/2001)>.
3. K. Uekama, S. Narisawa, F. Hirayama, M. Otagiri. *Int. J. Pharm.* **16**, 327 (1983).
4. I. C. R. Galán, R. López, A. C. Galán. *Ciência Técnica* **6**, 366 (1987).
5. K. Terada, K. Yamamoto, Y. Nakai. *Expo-Congr. Int. Technol. Pharm.*, 3rd, **5**, 246 (1983).
6. K. Uekama, F. Hirayama, A. Fujise, M. Otagiri, K. Inaba, H. Saito. *J. Pharm. Sci.* **73**, 382 (1984).
7. A. L. Thakkar, P. B. Kuehn, J. H. Perrin, W. L. Wilham. *J. Pharm. Sci.* **61**, 1841 (1972).
8. C. P. Rossel, J. S. Carreño, M. R. R-Baeza, J. B. Alderete. *Quím. Nova* **23**, 749 (2000).

Madalena Pinto
madalena@ff.up.pt

High standards in safety measures should be maintained in all work carried out in Medicinal Chemistry Laboratories. The handling of electrical instruments, heating elements, glass materials, dissolvents and other inflammable materials does not present a problem if the supervisor's instructions are carefully followed.

This document has been supervised by Prof. MADALENA PINTO (madalena@ff.up.pt) who has informed that no special risk (regarding toxicity, inflammability, explosions), outside of the standard risks pertaining to a Medicinal Chemistry laboratory exist when performing this exercise.

If your exercise involves any "special" risks, please inform the editor.