PRÁCTICA V.5
BENZODIAZEPINAS – PARTE I
AISLAMIENTO E IDENTIFICACIÓN EN FÁRMACOS COMERCIALES
BENZODIAZEPINES – PART I. ISOLATION AND IDENTIFICATION IN
COMMERCIAL DRUGS

Miriam A. Martins Alho, Mirta L. Fascio y Norma B. D'Accorso

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1. Analytical Techniques (Cap. Colour Test, H. M. Stevens; Cap. Thin-Layer Chromatography, A. C. Moffat).
2. Analytical Profiles of drug Substances, Ed. Klaus Florey. Academic Press, San Diego, CA USA, 1984
3. The Merck Index, An Encyclopedia of Chemicals, Drugs and Biologicals (1989)

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EXERCISE V.5
BENZODIAZEPINES – PART I
ISOLATION AND IDENTIFICATION IN COMMERCIAL DRUGS

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
E-mail: norma@qo.fcen.uba.ar

Introduction


The benzodiazepines are included in the “psychotropic agents” class. Some of these drugs (diazepam, clonazepam, oxazepam, bromazepam) are used for treatment of anxiety or insomnia disorders, muscle or other chronic pain (such as arthritis) in more or less proportion, as well as muscle relaxant and in epilepsy treatment. Other benzodiazepines such as nitrazepam or flurazepam are more selective drugs.

The basic structure is obtained by a condensation reaction between an *o*-aminobenzophenone derivative and an amino acid moiety; these compounds can be hydrolyzed to obtain the precursor drugs, from which we can infer the possible degradation mechanism of these kinds of medicines.


The most commonly used benzodiazepines in the medicinal field are:


1. Diazepam


2. Bromazepam


3. Flunitrazepam





4. Flurazepam

5. Lorazepam

6. Nitrazepam

**7. Oxazepam****8. Prazepam****9. Nordazepam****10. Clonazepam**

The objective of this experiment is the isolation and identification of the active compound from a commercial medication. The different methodologies we propose depend on the equipment facilities in each laboratory. The experiment is divided into three parts.

Experimental Procedure

1. Extraction

Crush four or five tablets of the unknown commercial medication in a mortar and extract the resulting powder with methylene chloride (2 x 15 mL). Filter the mixture through a celite bed to remove the insoluble compounds (mainly excipients). The obtained colorless solution will be used for the assays.

2. Thin-layer and preparative chromatography

As it is known (of course), the R_f chromatography value cannot be used as an identification criterion, but in a series or families of compounds, such as pesticides and alkaloids, among others, the use of standard chromatographic solvents allows us to compare with standard data. Particularly, for benzodiazepines, the solvent systems are:

Ethyl acetate [TF]

Chloroform:acetone (4:1) [TD]

Ethyl acetate:methanol:ammonia (85:10:5) [TE]

Perform thin-layer chromatography of the previously obtained extract (item 1) using three chromatographic systems mentioned above and compare the R_f values with the bibliographic data listed below.

Table 1

BENZODIAZEPINE	[TF]	[TD]	[TE]
1. Diazepam	0.22	0.15	0.44
2. Bromazepam	0.20	0.13	0.62
3. Flunitrazepam	0.48	0.58	0.77
4. Flurazepam	0.03	0.03	0.74
5. Lorazepam	0.39	0.23	0.45
6. Nitrazepam	0.45	0.35	0.59
7. Oxazepam	0.37	0.22	0.44
8. Prazepam	0.55	0.64	0.81
9. Nordazepam	0.45	0.34	0.67
10. Clonazepam	0.45	0.35	0.60

Because the R_f is not an identification criterion, it is necessary to isolate the active drug for confirmation, particularly when the R_f values are similar. If there is only one stable drug, the active compound can be isolated directly from the item 1 extract. Otherwise, the mixture of components must be separated by preparative thin-layer chromatography, using ethyl acetate:acetone (4:1). After the chromatographic plate is developed, the spots are located by exposing the border of the slide to iodine vapor. The portion of silica identified is removed, and the benzodiazepines are extracted with dichloromethane, avoiding the strip that was in contact with the iodine. The extract is filtered through a celite bed, and the solvent is evaporated under diminished pressure. The resulting solid is dried, and its melting point is measured. The benzodiazepine identity is assured by means of a mixed melting point using appropriated standards.

Table 2

BENZODIAZEPINE	Fusion point (°C)
1. Diazepam	125–126
2. Bromazepam	237–238.5
3. Flunitrazepam	166–167
4. Flurazepam	77–82
5. Lorazepam	166–168
6. Nitrazepam	224–226
7. Oxazepam	205–206
8. Prazepam	145–146
9. Nordazepam	216–217
10. Clonazepam	236.5–238.5

3. Functional group characterization

When the correct identity of the benzodiazepine cannot be established or does not have the standard to make the mixed melting point, a complementary test must be performed.

According to the equipment facilities, two alternative techniques can be used:

- a. EI-MS
- b. Functional group analysis

a. All commercial benzodiazepines can be identified by this method because their

molecular weights are smaller than 400 amu. A simple mass spectrometer can be used to detect the molecular ion. When bromine or chlorine are present, the isotopic clusters can be observed for all halogenated fragments.

b. This option is suitable to distinguish between most of the commercial drugs. Exception must be made for lorazepam vs. oxazepam (*i*) or prazepam vs. nordazepam (*ii*), but the components of both couples have great differences between their melting points. Despite that, these assays can be made for confirmation:

i. A quantitative chlorine gravimetry to distinguish lorazepam from oxazepam, which have similar R_f .

ii. If any doubts remain in distinguishing prazepam from nordazepam by R_f values in different solvent systems (**TD** or **TE**), FTIR can be performed. Differentiation can be made by looking for a band corresponding to N-H stretching. Another option is to perform a ^1H NMR spectrum and identify the signal corresponding to the cyclopropylmethyl group of prazepam.

Assays for hydroxyl and nitro groups can be made directly over the sample, but also by using FTIR, which is the safest way to see these functional groups. To analyze the presence of halogens in the molecule without a mass spectrometer, a sample mineralization must be made, using sodium fusion.

Sodium fusion 1

With a pair of tweezers, remove a sodium metal piece and dry it carefully with dry filter paper. Cut a little cube of sodium (not more than 1–2 mm by side) and place it in a clean and dry 10 × 75 mm Pyrex test tube. Heat the tube with a Bunsen burner. When about 1 cm of the tube is filled with sodium vapors, put it on vertical and add one drop of liquid or 5 mg of unknown solid. This addition can produce flame and smoke. Heat the tube for another 60 seconds and allow it to cool to room temperature. Place the tube in a 50 mL beaker and add 1 mL of alcohol to destroy the sodium excess, aided with a glass rod. If the tube breaks, the ethanolic solution can be recovered from the beaker. Add 2 mL of water and boil the mixture for a few seconds, and filter the colorless mixture. Cool this clear filtrate and use it in the following tests.

Note: if the filtrate is not colorless, the quantity of sample was too much, and the fusion must be repeated.

CAUTION: Metallic sodium must be manipulated with great care. Never use large quantities of metallic sodium, and do not touch the sodium with your fingers, use tongs or pincers. Do not throw small residual pieces of sodium into water, or into a sink or waste crock; destroy sodium residues in an appropriate quantity of alcohol.

Carry out the sodium fusion with care, and be especially cautious in decomposing the fused mass with alcohol. If too large an amount of sodium is used, the resulting treatment with alcohol may be too vigorous. Safety glasses must be worn for the experiment and should be worn at all times in the laboratory.

General halogen test

Use 6 M nitric acid to acidify 1 mL of the stock solution. Boil the acidified aliquot gently for 5 min in the hood to remove the hydrogen cyanide formed during the mineralization. One or two drops of silver nitrate solution are added to 0.5 mL of solution. Precipitate formation indicates the presence of halogens, except fluoride. The solid is allowed to precipitate and separate from the solution. 0.5 mL of diluted ammonium hydroxide (50 % in volume) is added to the solid; if dissolution is observed, the halide is chloride, if there is no dissolution, it is bromide.

To determine fluoride: acidify 1 mL of the fusion solution with acetic acid, boil, and then when the solution arrives at room temperature, a filter paper moistened with zirconium-alizarin solution [3 mL of alizarin solution in ethanol (1 %) to 2 mL of zirconium chloride or zirconium nitrate (0.4 %) test the evolved gases] is introduced. The paper is submerged and then dried. A drop of acetic acid (50 %) is added to humidify, before use. Put a drop of acidified solution on the treated paper, a yellow coloration indicates the presence of fluoride.

Aromatic nitro group's characterization

The benzodiazepines are unstable in acid medium, due to a hydrolysis reaction. The identification of the nitro groups is based on the transformation of this group in the corresponding aniline, so the hydrolysis does not interfere with this determination. 10 mg of the substance is dissolved in 0.5 mL of acetic acid (heat if necessary). 0.25 mL of concentrate chloride acid and 0.6 mL of water are added, and then 0.2 g of Zn in powder is added carefully. The mixture is heated at reflux during 5 min, filtered, and the solvent is evaporated. Then a drop of β -naphthol solution (2 % in ethanol) is added. The red or orange indicates that nitro groups are present.

Hydroxyl group's determination

Five to ten drops of saturated solution of KMnO_4 in 4 N H_2SO_4 are added to a little portion of substance in the study, and are shaken for 2–3 min. If the MnO_2 precipitation appears in 5–7 min, this indicates that the benzodiazepines have hydroxyl groups (if the coloration takes a longer time to appear, this can be due to interferences present in the substances).

Test	1	2	3	4	5	6	7	8	9	10
Zr/aliz.	-	-	+	+	-	-	-	-	-	-
Cl^-	+	-	-	+	+	-	+	+	+	+
Br^-	-	+	-	-	-	-	-	-	-	-
NO_2^-	-	-	+	-	-	+	-	-	-	+
HO^-	-	-	-	-	+	-	+	-	-	-

References

1. Analytical Techniques (Cap. Colour Test, H. M. Stevens; Cap. Thin-Layer Chromatography, A. C. Moffat).
2. Analytical Profiles of Drug Substances, Ed. Klaus Florey. Academic Press, San Diego, CA USA, 1984
3. The Merck Index, An Encyclopedia of Chemicals, Drugs and Biological (1989).

Miriam A. Martins Alho, Mirta L. Fascio, and Norma B. D'Accorso

norma@go.fcen.uba.ar

High standards in safety measures should be maintained in all work carried out in Medicinal Chemistry Laboratories.

The handling of electrical instruments, heating elements, glass materials, dissolvents and other inflammable materials does not present a problem if the supervisor's instructions are carefully followed.

This document has been supervised by Prof. Dra. Norma B. D'Accorso (norma@go.fcen.uba.ar) who has informed that no special risk (regarding toxicity, inflammability, explosions), outside of the standard risks pertaining to a Medicinal Chemistry laboratory exist when performing this exercise.

If your exercise involves any "special" risks, please inform the editor.