**PRÁCTICA VI.6**  
**ANÁLISIS DE FLAVONOIDES EN PLANTAS**  
**ANALYSIS OF FLAVONOIDS IN PLANTS**

**Olga Lock, Isabel Cabello, Víctor Hugo Doroteo**

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
\*E-mail: olock@pucp.edu.pe

#### 1. INTRODUCTION

Flavonoids belong to a group of natural compounds that contain 15 carbon atoms in their basic nucleus and arranged in a C<sub>6</sub>-C<sub>3</sub>-C<sub>6</sub> configuration, that is, two aromatic rings linked by a three-carbon unit which may or may not be part

of a third ring; conventionally, these rings are labelled as A, B, and C.

These compounds occur as aglicones or O-glycosides in which one or more of their hydroxyl groups is bonded to one or more sugar molecules by an acid-labile hemiacetal bond. Sugars such as glucose, galactose, rhamnose, xylose, and arabinose are the most commonly involved, having a much higher possibility of being in C<sub>3</sub> or C<sub>7</sub> position.

Although flavonoids have been mainly used as wool colorants, some of their actions in plants involve (a) protecting them against UV and visible light exposure; (b) attracting bees, butterflies, and other insects in the pollination process; (c) repelling insect attacks; (d) and their behavior as growth regulators and enzyme inhibitors in vitro. Some flavonoids also exhibit a variety

of biological effects such as antiinflammatory, antiallergenic, antiulcerogenic, antiviral, anticarcinogenic, or antioxidant; besides that, they can be used in the treatment of diseases such as diabetes, cardiac problems, or in capillary fragility, among others.

Flavonoids possessing in their structure more than one unsubstituted hydroxyl groups, or a sugar, are polar compounds and because of this, are soluble in solvents such as ethanol, methanol, or water. In the same way, less polar aglycones tend to be more soluble in solvents such as ether and chloroform. Both groups usually show intense absorptions in the UV–vis region of the electromagnetic spectrum due to the presence of an aromatic conjugated system in their structures.

A preliminary classification of the flavonoid type in a plant extract can be performed, based on its solubility and color reaction behavior, followed by a chromatographic test.

In this experiment, we describe the tests for the detection of flavonoids in a plant extract using an easy color reaction (Shinoda reaction) and a thin-layer chromatographic analysis with suitable spray reagents, and the quantitative analysis by UV–vis spectroscopy technique with shift reagents.

## 2. EXPERIMENTAL PROCEDURE

### 2.1. QUALITATIVE ANALYSIS

#### 2.1.2. Sample preparation

Weigh 1 g of dried and powdered sample in a 50-mL Erlenmeyer flask, add 5 mL of methanol, and heat the mixture in a water bath for about 10 min at 60 °C (use a sonicator if possible). Filter the solution through a filter paper (A extract).

### 2.1.2. Shinoda reaction

Place in a test tube 20 drops of the A extract, add a few magnesium shavings, and some drops of concentrated hydrochloric acid. Look at the color change.

### 2.1.3. Thin-layer chromatography (TLC) analysis

- a. Apply 30  $\mu\text{L}$  of the A extract in a  $10 \times 3$  cm silica gel 60F<sub>254</sub> (Merck) TLC plate.
- b. Let the chromatography chamber saturate (use a beaker of 8 cm diameter and 11 cm height) with a mixture of ethyl acetate:methanol:water (100:13.5:10). Alternatively, a mixture of ethyl acetate:acetic acid:formic acid:water (100:11:11:26) can be used as mobile phase
- c. Place the plate in the chamber and let it develop until the mobile phase moves up to 8 cm height. Dry the plate with hot air (by means of a hair dryer) and detect the spots with long-wave UV–vis lamp and/or spraying the plate with a NP-PEG solution. Observe the fluorescence colors under UV light after spraying it.
- d. Calculate the  $R_f$  value for the marked spots.

## 2.2. QUANTITATIVE ANALYSIS BY UV–VIS SPECTROSCOPY

### 2.2.1. Sample preparation

- a. Weigh 1 g of dried and powdered sample in a 50-mL Erlenmeyer flask, add 20 mL of 80 % ethanolic aqueous solution. Place it in a sonicator for about 30 min at 60 °C (use a hot plate). Filter the solution into a 25-mL volumetric flask and bring up to volume with 80 % ethanolic aqueous solution.
- b. Place 100  $\mu\text{L}$  of the extract solution obtained in (a), in a 10-mL volumetric flask, add 200  $\mu\text{L}$  of 1 M potassium acetate and 200  $\mu\text{L}$  of 10 % aluminum nitrate solutions, bring up to volume with 80 % ethanolic aqueous solution.

- c. Let the solution stand for 40 min and determine the absorbance A at 415 nm.

### 2.2.2. Standard calibration curve

- Weigh 2.7 mg of quercetin in a 10-mL volumetric flask and dissolve with 80 % ethanolic aqueous solution (stock solution).
- Prepare four solutions, by placing 700.0, 350.0, 175.0, and 100.0  $\mu\text{L}$  from the stock solution in 10-mL volumetric flasks and adding to each one 200  $\mu\text{L}$  1M potassium acetate and 200  $\mu\text{L}$  10 % aluminum nitrate solutions. Bring up to volume with 80 % ethanolic aqueous solution.
- Determine A at 415 nm after 40 min for each solution and draw the standard curve.

### 2.2.3. Results

- Calculate the total flavonoid concentration using the standard curve, and express it as quercetin.
- As reference information, we are including the values for celery and chamomile. These contents can vary as a function of the sample source.

|                         | Celery | Chamomile |
|-------------------------|--------|-----------|
| Total flavonoids (mg/g) | 0.041  | 5.57      |

## 3. Notes

- The qualitative analysis is very simple to be followed. It is suitable to use aluminum foil plates which are, commercially available and can be cut to the desired size, for TLC analysis.

- b. The solvents used in the experiment are the ones commonly used in the laboratory; as an alternative spray reagent, a ferric chloride ethanolic solution may be used.
- c. The chromogenic agent NP-PEG is constituted using the following solutions: (a) 1 % diphenylboryloxyethylamine (NP) methanolic solution and (b) 5 % polyethylenglycol-400 (PEG) ethanolic solution. To use it, first spray (a), then (b).
- d. We consider the UV–vis spectroscopy method as the most suitable for the flavonoid quantification, because of its accessibility and availability of reagents.

## REFERENCES

1. Lock, O. 1994. *Investigación Fitoquímica, Método en el Estudio de Productos Naturales*. Fondo Editorial PUCP. Lima.
2. Mabry, T.J., Markham, K.R. and Thomas, M.B., 1970. *The Systematic Identification of Flavonoids*. Springer–Verlag, Berlín.
3. Oliveira, C.M., et. al. 1999. *Farmacognosia, da Planta ao Medicamento*. Editora DAUFSC, Florianapolis.
4. Wagner, H., Bladt, S., Zgainski, E.M., 1996. *Plant Drug Analysis. A Thin Layer Chromatography*. Springer-Verlag, Berlín.
5. Park, Y., Koo, M., Ikegaki, M., Contado, J., 1997. Comparison of the flavonoid aglycone contents of *Apis mellifera* propolis from various regions of Brazil. *Arq. Biol. Tecnol.* **40** (1), 97-106.

**Olga Lock**

olock@pucp.edu.pe

High standards in safety measures should be maintained in all work carried out in Medicinal Chemistry Laboratories.

The handling of electrical instruments, heating elements, glass materials, dissolvents and other inflammable materials does not present a problem if the supervisor's instructions are carefully followed.

This document has been supervised by Prof. OLGA LOCK SING ([olock@pucp.edu.pe](mailto:olock@pucp.edu.pe)) who has informed that no special risk (regarding toxicity, inflammability, explosions), outside of the standard risks pertaining to a Medicinal Chemistry laboratory exist when performing this exercise.

If your exercise involves any "special" risks, please inform the editor.